

Development of an In Vitro Diagnostic Multivariate Index Assay (IVDMIA) for Aggressive Prostate Cancer

EDRN Biomarker Reference Laboratory (BRL) at Johns Hopkins



Daniel W. Chan, Ph.D., DABCC (PI)

Lori Sokoll, Ph.D., Zhen Zhang, Ph.D., Hui Zhang, Ph.D.

June 30, 2020

JHU BRL Aims

The goals of the JHU BRL:

- **Collaborative studies:** Help accelerate the translation of biomarkers into clinical diagnostics by providing resources and support for the analytical and clinical validation of biomarker.
- **Product Development:** Develop an IVD MIA (*In vitro diagnostic multivariate index assay*) that combines a panel of serum biomarkers into a single-valued numerical index score. It is intended to evaluate the risk of aggressive prostate cancer in men prior to prostate biopsy.

JHU BRL Milestones Prostate Cancer


Candidate Biomarkers – Prostate					
Candidate Biomarker	Discovery			Pre-validation	Validation
	Discovery	Predictive Analysis	Assay Refinement	Blinded Limited Cross-Sectional	Large Cross-Sectional
BRL-IVDMIA for aggressive prostate cancer					
BRL-MiCheck, Minomic Ltd. Collaboration CLIA/LDT					
<i>Collaborative projects- see reports from EDRN lab PIs</i>					
Sanda CVC – Biomarker Panel – phi and PCA3					
Sanda CVC – phi in PASS study					
UTHSC Liss CVC – phi and imaging					
Collaborative group study – Imaging and Biomarkers (PI Wei)					

JHU BRL Milestones

Ovarian, Pancreas, and Liver Cancers


	Discovery			Pre-validation	Validation
	Discovery	Predictive Analysis	Assay Refinement	Blinded Limited Cross-Sectional	Large Cross-Sectional
Candidate Biomarkers - Ovarian					
<i>Collaborative projects- see reports from EDRN lab PIs</i>					
Collaborative group study – Skates, Early detection biomarker panel (PI Skates)					
Collaborative group study – Uterine lavage biomarker panel study (PI Skates)					
Candidate Biomarkers - Pancreas					
BRL-Multiplex Immunoassay panel					
Candidate Biomarkers - Liver					
BRL - Biomarkers panel, Abbott Labs collaboration					

JHU BRL Product Development Biomarker Identification (1)




JHU BRL Product Development Biomarker Identification (2)

TIE-2 (soluble angiopoietin-2)




Danni Li, Hanching Chiu, Vinita Gupta, Daniel W. Chan.
Clinica Chimica Acta, (2012)413:1506–1511

Fucosylated PSA




Qing Kay Li, Li Chen, MH Ao, Joyce Chiu, Zhen Zhang, Hui Zhang, Daniel W Chan.
Theranostics (2015) 5:267-276.




Wang C, Hôti N, Lih TM, Sokoll LJ, Zhang R, Zhang Z, Zhang H, Chan DW.
Clinical Proteomics. (2019) Apr 6;16:13.


JHU BRL Product Development Biomarker Identification (3)




JHU BRL Product Development – Pre-validation




JHU BRL Product Development Biomarker Distribution



Panel Selection (using backward stepwise multivariate logistic regression): complement *phi* to detect AG from NAG Pca (from bootstrap estimation)




Panel Selection (using backward stepwise multivariate logistic regression): complement phi to detect AG from NAG PCa + non PCa (from bootstrap estimation)



Improvement in specificity at 95% sensitivity

	AUC (95% CI)	SN (%)	SP (%)	Neg	Pos	False-Neg	False-Pos
AG vs NAG							
Panel 1	0.942 (0.876-1.000)	95.0	76.0	19	24	1	6
<i>phi</i> + Fuc-PSA	0.914 (0.828-0.980)	95.0	76.0	19	24	1	6
<i>phi</i>	0.872 (0.748-0.971)	95.0	56.0	14	24	1	11
PSA	0.866 (0.749-0.956)	95.0	44.0	11	24	1	14
AG vs NAG & Non							
Panel 2	0.934 (0.866-0.987)	95.0	78.2	43	24	1	12
<i>phi</i> + Fuc-PSA	0.918 (0.842-0.974)	95.0	69.1	38	24	1	17
<i>phi</i>	0.898 (0.814-0.963)	95.0	65.5	36	24	1	19
PSA	0.807 (0.697-0.905)	95.0	36.4	20	24	1	35

Summary - JHU BRL Product Development



Acknowledgments

JHU BRL Team

Daniel W. Chan, Ph.D.

Lori Sokoll, Ph.D.

Zhen Zhang, Ph.D.

Hui Zhang, Ph.D.

Jin Song, MD, Ph.D.

Shiyong Ma, Ph.D.

Naseruddin Höti, Ph.D.

Rodrigo V. Eiguez

Phaedre Mohr

Renu Dua

Emory CVC Team

Martin G. Sanda, MD

Datta Patil

Kristen Douglas May

Sierra Williams

Rebecca Arnold

NCI

Jacob Kagan, Ph.D.

Lynn Sorbara, Ph.D.

Development and evaluation of the MiCheck test for aggressive prostate cancer

- MiCheck is an IVD/MIA for aggressive prostate cancer using a logistic regression model with the following variables: PSA, patient age, DRE, Leptin, IL-7, VEGF, Glypican-1.
- *Shore ND, Pieczonka CM, Henderson RJ, Bailen JL, Saltzstein DR, Concepcion RS, Beebe-Dimmer JL, Ruterbusch JJ, Levin RA, Wissmueller S, Le TH, Gillatt D, Chan DW, Campbell DH, Walsh BJ. Urol Oncol. 2020 Apr 16:S1078-1439(20)30097-1.*
- Developed by Minomic Ltd.
- Agreement with Cirrus Dx, a CLIA lab, MiCheck performed as an LDT (lab developed test).

The MiCheck® Prospective Trial


- Principal Investigator: Dr Neal Shore
- Prospective, non-randomized case-control study at 12 US research centers, 320 subjects
- All subjects proceeding to Prostate Biopsy on the basis of an elevated, age-adjusted PSA score
- Total 320 patients falling in to three groups as follows:
 - Group A: 141 patients (44%), no cancer
 - Group B: 62 patients (19%), GS = 3+3, CaP
 - Group C: 117 patients (37%), GS \geq 3+4, CaP
- Serum + plasma samples collected
- Standardized PSA test + centralized pathology review by Dr Scott Lucia (both Gleason and ISUP scores)
- Algorithm developed for differentiation of aggressive (GS \geq 3+4) v non-aggressive cancer (GS=3+3)/ no cancer



Past President
Neal D. Shore, MD, FACS
Myrtle Beach, SC

Study designed by Minomic and its International clinical advisory panel

Key MiCheck® Prostate Test Outcomes



- ✓ Highest AUC when directly compared to PSA, PHI and % free PSA
 - ✓ Best test specificity at an extremely high test sensitivity (95%)
 - ✓ Results in 43% fewer unnecessary biopsies
 - ✓ High Negative Predictive Value (97%)
- MiCheck® has demonstrated ability to differentiate aggressive from non-aggressive cancer and no cancer
 - MiCheck® demonstrates a compelling advantage in specificity over existing products shown here
 - The other test's lack of specificity has led to criticism due to the resulting unnecessary biopsies – clinicians uniformly desire high specificity
 - The MiCheck® technology helps fill this substantial unmet clinical need by allowing a reduction in number of biopsies

When finding aggressive cancer MiCheck® wins hands down